

From the makers of MyClyns®  
Personal Protection Spray

**Clyns**<sup>™</sup>  
BRANDS

Germ Defense Solutions

# ViralClyns<sup>™</sup>

***ViralClyns<sup>™</sup> is the only cost-effective, flat-fold N-95 respirator with the comfort and breathability of a face mask plus an infused EPA-approved antimicrobial.***

***N-95 Disposable Respirator***  
***with EPA-registered Antimicrobial Protection***


- *NIOSH-certified N95 respirator*
- *Easy to use, simple to store*
- *Effective against all types of microorganisms - including bacteria, viruses, fungi, bacterial spores, and protozoa*
- *Most protective and cost-effective N95*
- *Available in two sizes*



# The ViralClyns Difference

	<b>ViralClyns</b> N95 Respirator	<b>N95</b> Respirator	<b>Surgical</b> Mask
<b>EPA-Approved Antimicrobial Layer</b>	✓		
<b>N-95 Rated</b>	✓	✓	
<b>Breathability</b>	✓		✓
<b>Easy to Store</b>	✓		✓
<b>Comfort</b>	✓		✓
<b>Face Fit</b>	✓	✓	
<b>Low Inward Leakage</b>	✓		

## Engineered for Better Protection



From its antimicrobial-coated scrim and protective filtration layer to its secure face seal and elastic head straps, the ViralClyns N95 Respirator with Antimicrobial Protection is designed to meet your needs. There are no compromises when it comes to comfort, face fit, breathability, and the added feature of an antimicrobial. You won't find another respirator like this one.

Contact us at ▶▶▶▶

# 1-877-Go-Clyns

or for more information ▶▶▶▶

## www.ClynsBrands.com

† NIOSH has certified ViralClyns respirators to meet the requirements in 42 CFR 84.

## When is the ViralClyns Used?

- Crime Scene
- Accident Scene
- Cold/Flu Season
- Pandemics/Epidemics
- Body Pick-up
- Disaster Response
- Biological Threats

### AVAILABLE IN TWO SIZES

**S/M**

**M/L**

### A430/A450 NIOSH Test Results Sodium Chloride (NaCl) N95 Filtration

# Tested	Flow Rate (L/min)	Average Initial Filter Resistance (mm of H <sub>2</sub> O)	Maximum Allowable % Leakage	Initial % Leakage Average (± SD)	Maximum % Leakage Average (± SD)	Result Pass or Fail
20	85	9.0	5	1.142 (± 0.3819)	1.247 (± 0.4136)	PASS

### A430/A450 NIOSH Test Results Exhalation and Inhalation Resistance

	# Tested	Maximum Allowable Resistance (mm of H <sub>2</sub> O)	Actual Resistance (mm of H <sub>2</sub> O) Average (± SD)	Result Pass or Fail
Exhalation Resistance	6	25	5.3 (± 0.16)	PASS
Inhalation Resistance	6	35	5.3 (± 0.24)	PASS